

Apricot Personal Pipettors™ (PP) Multistation Selection Guide

	PP-96H-M	PP-96L-M	PP-96H-H	PP-96S-100	PP-384-M
No. of Labware Stations	6	6	6	6	6
No. of Pipetting Channels	96	96	96	96	384
Max Capacity per Channel (ul)	550	125	550	100	125
96-well	√	√	√	√	√
384-well	√	√	√	√	√
1536-well					√
Serial Dilution col to col	√	√	√		√
Serial Dilution row to row	√	√	√		√
SPE			√		

Designs that make sense

Apricot Designs, Inc. 681 Arrow Grand Circle, Covina, CA 91722, USA Tel: (626) 966-3290

