

Apricot Total Pipetting Solution (TPS) Selection Guide

	TPS-24	TPS-96H	TPS-96L	TPS-384
No. of Labware Stations	6	12,15	12,15	12,15
No. of Pipetting Channels	24	96	96	384
Max Capacity per Channel (ul)	125	550	125	125
96-well	√	√	√	√
384-well	√	√	√	√
1536-well	√			√
Serial Dilution col to col	√	√	√	√
Serial Dilution row to row	√	√	√	√

Designs that make sense

Apricot Designs, Inc. 681 Arrow Grand Circle, Covina, CA 91722, USA Tel: (626) 966-3290

